

Gastric Lavage & Tubing



Ethox Medical Products Division
International

www.ethoxint.com

Tum-E-Vac® Gastric Lavage Kits

Ethox International's Tum-E-Vac® is a pre-assembled, completely closed system designed to protect healthcare practitioners from possible exposure to stomach contents and charcoal solutions during gastric lavage or decontamination procedures. The Tum-E-Vac is a single-use, disposable kit that promotes quick, efficient lavage of stomach contents.

Product Features

Safe

- Disposable, single patient use
- Closed system prevents exposure to stomach contents and charcoal solutions
- Latex free

Versatile

- Available with Activated Charcoal suspension with or without Sorbitol
- Twist Lok™ expressing key to aid in charcoal administration
- Easy fill top port
- Available with Universal Solution Spike for use with pre-filled solution containers
- Pre-assembled, quick and easy to use
- Kits contain 3L solution bag and 5L drain bag



Gastric Lavage Kits without Charcoal

Catalog No.	Product	Quantity
2055	Lavage Kit with 3L solution bag	10 units per box
2060	Bite Block	20 units per box
2075	Lavage Kit with spike, no solution bag	10 units per box
3010	Actidoes-Aqua - 25g tube	12 units per box
3020	Actidoes-Sorbitol - 25g tube	12 units per box

ISO Certified
CE Marked



Lavage Convenience Kits with Lavage Tubes and/or Charcoal

Catalog No.	Product	Quantity
2054	Kit w/ 32 fr Lavage Tubes, 3L bag	10 units per box
2057	Kit w/ 32 fr Lavage Tubes, (2) 25g Activated Charcoal	10 units per box
2058	Kit w/ 32 fr Lavage Tubes, (2) 25g Activated Charcoal w/ Sorbitol	10 units per box
2059	Kit w/ 32 fr Lavage Tubes. (1) 25g Act. Charcoal, (1) 25g Act. Charcoal w/ Sorbitol	10 units per box
2078	Kit w/ spike, 32 fr Lavage Tubes, (2) 25g Activated Charcoal w/ Sorbitol	10 units per box
2081	Kit w/ spike, (2) 25g Activated Charcoal w/ Sorbitol	10 units per box

Lavage Tubes

Ethox International's Lavage Tubes have an open bullet tip design with multiple eyelets to create free-flowing movement with low resistance. All tubes are individually packaged, clean and ready to use.

Product Features

- Disposable, single patient use
- Soft clear tube with smooth tip
- Suction connector on proximal end
- Radiopaque stripe to help confirm proper tube placement
- Latex free

Lavage Tubes

Latex free, disposable, radiopaque, open bullet tipped design, multiple eyes. 48" long with reference marks at 18", 22", and 20" from the distal end.

Catalog No.	Product	Quantity
2124	24 fr Lavage Tube	10 units per box
2128	28 fr Lavage Tube	10 units per box
2132	32 fr Lavage Tube	10 units per box
2136	36 fr Lavage Tube	10 units per box



500-709 0810