

# Wound Irrigation



**Ethox** Medical Products Division  
International

www.ethoxint.com


## Combiguard II®

### Irrigation Splash Guard

Ethox International's Combiguard II® is a sterile, economical solution to prevent wound irrigation splashback. The Combiguard II produces a fluid stream equivalent to that of an 18g needle and attaches to luer tip syringes. The clear durable splash guard provides protection for the practitioner while the dual outlet ports provide an escape path for the contaminated fluids.

#### Product Features

- Egg-shaped design ideal for larger coverage area
- Rounded edge provides smooth surface against the skin
- Angled fluid delivery inlet directs fluid away from cleansed area and toward outlet ports
- Dual outlet ports provide escape path



## Combiport®/Combiguard® Kit

### Wound Irrigation Device

Combiport® is an effective instrument for administering a high volume of fluid in a short period of time. Patented spike attaches to standard 500ml or 1000ml normal saline bag. Irrigation removes foreign bodies, reduces pathogen numbers, and moistens the wound's tissue.

#### Product Features

- Convenient sterile peel open pouch
- Removable cap for easy stop/start during cleansing
- Patented end port fits standard IV bags
- Cost effective with minimal set up time



## Combiguard II & Combiport

Catalog No.	Product	Quantity
10100	Combiport®/Combiguard® Kit (includes Combiguard®)	50 per case
10222	Combiguard II® for use with 30-40cc luer tipped syringe	50 per case
10222NS	Combiguard II® - Non-Sterile	500 per case

**ISO Certified  
CE Marked**