

Surgi-Stuf® Surgi-Kit® Lap-Kit®



Surgi-Stuf® Hand/Foot Pack


The polyester fiber that is used to make Ethox International's Surgi-Stuf® is easy to tear and aids in the precise positioning of the feet or hands by separating fingers and toes for a variety of medical procedures. The product is made from a voluminous material that allows pressure to diffuse from the hands or from the feet preventing potential tissue constriction. The pack is nonabsorbent which means that it remains supportive and will provide sufficient compression. Surgi-Stuf® is hypo-allergenic and sterile.



Surgi-Kit® Instrument Holders and Suture Organizers

These disposable instrument holders are available in multiple configurations to accommodate surgeon or physician needs. With the use of the Surgi-Kit®, instruments can be kept close to the operating field, which in turn saves time. The kit attaches securely to the operating drape with either adhesive tape or surgical clamps and the instruments are securely held in place with cords or hoses. The kits are sterilized, therefore eliminating the need for cleaning procedures.

*Also Available in Non-Sterile Bulk



Lap-Kit® Laparoscopic Instrument Holders

Like the Surgi-Kit®, Ethox International's Lap-Kit® offers the same convenience and quality.

Lap-Kit® has an extended length to accommodate all instruments up to 30 centimeters long and is useful for laparoscopic surgeries and laser procedures.

Order Information



Surgi-Stuf® Hand/Foot Pack

8001 Hand/Foot Pack - 12"L x 8.5"W - 25 units per Box



Surgi-Kit® Instrument Holders and Suture Organizers

- 100-SK 8 Pkt without tape-Sterile - 19cm x 23cm - 50 per Box
- 100-SK 8 Pkt without tape -Sterile - 19cm x 23cm - Case of 200
- 100N-SK 8 Pkt without tape -NonSterile - 19cm x 23cm - Case of 500
- 150-SK 4 Pkt without tape -Sterile - 19cm x 11cm - 100 per Box
- 150-SK 4 Pkt without tape -Sterile - 19cm x 11cm - Case of 400
- 150N-SK 4 Pkt without tape -NonSterile - 19cm x 11cm - Case of 500
- 200-SK 2 Pkt without tape -Sterile - 28cm x 18cm - 50 per Box
- 200-SK 2 Pkt without tape -Sterile - 28cm x 18cm - Case of 200
- 200N-SK 2 Pkt without tape -NonSterile - 28cm x 18cm - Case of 500
- 300-SK 8 Pkt with tape-Sterile - 19cm x 23cm - 50 per Box
- 300-SK 8 Pkt with tape-Sterile - 19cm x 23cm - Case of 200
- 300N-SK 8 Pkt with tape-Non-Sterile - 19cm x 23cm - Case of 500
- 350-SK 4 Pkt with tape-Sterile - 19cm x 11cm - 100 per Box
- z350-SK 4 Pkt with tape-Sterile - 19cm x 11cm - Case of 400
- 350N-SK 4 Pkt with tape-Non-Sterile - 19cm x 11cm - Case of 500
- 400-SK 2 Pkt with tape-Sterile - 28cm x 18cm - 50 per Box
- 400-SK 2 Pkt with tape-Sterile - 28cm x 18cm - Case of 200
- 400N-SK 2 Pkt with tape-Non-Sterile - 28cm x 18cm - Case of 500
- 650-SK 10 Pkt with tape-Sterile - 15cm x 19cm - 50 per Box
- 650-SK 10 Pkt with tape-Sterile - 15cm x 19cm - Case of 200
- 650N-SK 10 Pkt with tape- Non-Sterile - 15cm x 19cm - Case of 500



Lap-Kit® Laparoscopic Instrument Holders

- 700-LK 2 Pkt with tape-Sterile - 68mm x 30cm - Box of 25
- 700-LK 2 Pkt with tape-Sterile - 68mm x 30cm - Case of 100
- 700N-LK 2 Pkt with tape-Non-Sterile - 68mm x 30cm - Case of 200
- 800-LK 4 Pkt with tape-Sterile - 33mm x 30cm - Box of 25
- 800-LK 4 Pkt with tape-Sterile - 33mm x 30cm - Case of 100
- 900-LK 3 Pkt with tape-Sterile - (1)68mm x 30 cm (2)33mm x 30 cm - Box of 25
- 900-LK 3 Pkt with tape-Sterile - (1)68mm x 30 cm (2)33mm x 30 cm - Case of 100
- 900N-LK 3 Pkt with tape-Non-Sterile - (1)68mm x 30 cm (2)33mm x 30 cm - Case of 200

