

# Surgi-Cuff®

## Adult and Neonatal Disposable Blood Pressure Cuffs



In an effort to reduce the transfer of infections, prevent harmful cross contamination, and to present a more efficient and economical means of blood pressure testing and observation – Ethox International offers the Surgi-Cuff® disposable blood pressure cuff.

These are individually packaged and cleaned cuffs that are ready for use. They are available in a full range of circumference sizes and models for varying connectors. The cuffs also utilize hook and loop fasteners to provide a snug and exact fit.

The Surgi-Cuff® meets the American Heart Association's recommendations for accuracy, therefore ensuring consistency and reliability.

The Surgi-Cuff® is also available for infants. The neonatal line of cuffs provides the same quality and cost-effective advantages.

The infant models are sterile and individually wrapped, preventing the introduction of infection to the undeveloped immune system of a neonate. The cuff is gentle and clear with a soft, latex free, vinyl design. This prevents skin irritations and allows the patient's skin to be visible during the procedure. The cuff is designed with soft rounded corners and edges for increased protection and comfort.



### Product Features

- **Versatile**
  - One universal style
  - Compatible with an extensive variety of electronic monitors and manual sphygmomanometers
  - Adapters are also available
  - Available in both single and double tube models
- **Convenient**
  - The soft vinyl surface wipes clean easily
  - Individually packaged
  - The cuffs are completely disposable
  - Solid case molded from glass fiber reinforced polypropylene
- **Safe**
  - No latex bladder
  - Meets the American Heart Association's Recommendations

Tel: 1.800.521.1022

**Ethox International, Inc.**  
ethoxint.com

Fax: 716.842.2546

# Order Information



## Surgi-Cuf® Adult Disposable Blood Pressure Cuffs

<b>Child Size - 15-22cm Arm Circumference</b>	
5106	1-Tube Female luer Lock Connector, 10 per box (2lbs)
5406	1-Tube Quick-disconnect Connector, 10 per box (2lbs)
5206	2-Tube Multi-purpose Screw Connector, 10 per box (2lbs)
<b>Small Adult Size - 17-25cm Arm Circumference</b>	
5107	1-Tube Female luer Lock Connector, 10 per box (2lbs)
5407	1-Tube Quick-disconnect Connector, 10 per box (2lbs)
5207	2-Tube Multi-purpose Screw Connector, 10 per box (2lbs)
<b>Adult Size - 24-32cm Arm Circumference</b>	
5108	1-Tube Female luer Lock Connector, 10 per box (2lbs)
5408	1-Tube Quick-disconnect Connector, 10 per box (2lbs)
5208	2-Tube Multi-purpose Screw Connector, 10 per box (2lbs)
<b>Adult Long Size - 28-37cm Arm Circumference</b>	
5108P	1-Tube Female luer Lock Connector, 10 per box (2lbs)
5408P	1-Tube Quick-disconnect Connector, 10 per box (2lbs)
5208P	2-Tube Multi-purpose Screw Connector, 10 per box (2lbs)
<b>Large Adult Size - 34-42cm Arm Circumference</b>	
5109	1-Tube Female luer Lock Connector, 10 per box (3lbs)
5409	1-Tube Quick-disconnect Connector, 10 per box (3lbs)
5209	2-Tube Multi-purpose Screw Connector, 10 per box (3lbs)
<b>Thigh Model - 42-50cm Arm Circumference</b>	
5110	1-Tube Female luer Lock Connector, 10 per box (4lbs)
5410	1-Tube Quick-disconnect Connector, 10 per box (4lbs)
5210	2-Tube Multi-purpose Screw Connector, 10 per box (4lbs)

## Surgi-Cuf® Neonatal Disposable Blood Pressure Cuffs

<b>Size 1 - 3-6cm Arm Circumference</b>	
5101	1-Tube Male Luer Slip Connector, Case of 10 (1lb)
5201	2-Tube Male Luer Slip, Case of 10 (1lb)
<b>Size 2 - 6-9.5cm Arm Circumference</b>	
5102	1-Tube Male Luer Slip Connector, Case of 10 (1lb)
5202	2-Tube Male Luer Slip, Case of 10 (1lb)
<b>Size 3 - 8-12cm Arm Circumference</b>	
5103	1-Tube Male Luer Slip Connector, Case of 10 (1lb)
5203	2-Tube Male Luer Slip, Case of 10 (1lb)
<b>Size 4 - 9-14cm Arm Circumference</b>	
5104	1-Tube Male Luer Slip Connector, Case of 10 (1lb)
5204	2-Tube Male Luer Slip, Case of 10 (1lb)
<b>Size 5 - 12-18cm Arm Circumference</b>	
5105	1-Tube Male Luer Slip Connector, Case of 10 (1lb)
5205	2-Tube Male Luer Slip, Case of 10 (1lb)



## Surgi-Cuf® Adult BP Cuff Adapters

5000B	Inflation Bulbs - 5 per Box
5002	1-Tube Adapters for HP Monitors, 10 per Box (1lb)
5004	2-Tube Adapter for Marquette Monitors, 10 per Box (1lb)
5006	2-Tube Adapters for Datex Monitors, 10 per Box

**Ethox International, Inc.**

Tel: 1.800.521.1022

ethoxint.com

Fax: 716.842.2546