

# Tum-E-Vac®

## Gastric Lavage



### Gastric Lavage Kit

Ethox International is proud to showcase the innovative Tum-E-Vac® gastric lavage kit. This a pre-assembled closed system that is designed to prevent exposure to stomach contents and charcoal solutions during gastric lavage or decontamination procedures.

This kit has all of the necessary components conveniently sealed and ready for use. The kit is easily utilized by a single person and creates an efficient means for evacuation. Additionally, the kit is completely disposable, which not only simplifies the traditional lavage process by eliminating cleaning procedures, but furthermore, it prevents contamination.

### Gastric Lavage Tubes

Ethox's lavage tubes are of superior quality and are designed for effective use. These lavage tubes have an open tip design with multiple eyelets to create free flowing movement with low resistance. This speeds up the procedure without having to use mucosa damaging suction. All tubes are individually packaged which prevents possibilities of contamination. There is no cleaning required and the product is ready for immediate use.

## Product Features

### Gastric Lavage Kit

- Available with Activated Charcoal Suspension in Specialized Soft Tubes
- Available with Universal Solution Spike
- Includes a disposable Bite Block and 600cc Piston Syringe
- Completely Disposable
- Multiple specialized configurations and options

### Lavage Tubes

- Latex free
- Created with a Radiopaque Stripe
  - Confirms proper tube placement
- Available in both single and double lumen models



## Order Information



### Tum-E-Vac®

#### Gastric Lavage Kits

Catalog Number	Description
2055	Lavage Kit w/ 3L bag - 10 units per Box (11lbs)
2045	Pediatric Lavage Kit - 10 units per Box (9lbs)
2075	Lavage Kit w/spike, no bag - 10 units per Box (10lbs)
2060	Bite Black - 20 units per Box (1lb)
3010	Actidoes-Aqua, 25g tube - 12 per Box
3020	Actidoes-Sorbitol, 25g tube - 12 per Box

### Tum-E-Vac®

#### Lavage Convenience Kits with Lavage Tubes and/or Charcoal

2054	32fr Lavage Tubes - Lavage Kit w/32fr & 3L bag - 10 units per Box
2057	32fr Lavage Tubes - Lavage Kit w/32fr & Aqua + 2-25g Activated Charcoal - 10 per Box (19lbs)
2058	32fr Lavage Tubes - Lavage Kit w/32fr & Sorbital + 2-25g Activated Charcoal w/Sorbitol - 10 per Box (19lbs)
2059	32fr Tubes - Lavage Kit w/32fr, Sorbital & Aqua + 1-25g Activated Charcoal & 1-25g Activated Charcoal w/Sorbitol - 10 per Box (19lbs)
2078	32fr Lavage Tubes - Lavage Kit w/spike & 2 Sorbitol + 2-25g Activated Charcoal w/Sorbitol - 10 per Box (19lbs)
2081	Lavage Kit w/ Spike, Sorbitol & Aqua + 2-25g Activated Charcoal w/Sorbitol - 10 per Box (12lbs)



### Bite Blocks

Will accommodate up to a 40 fr. tube. Latex free, individually packaged, and completely disposable.

2060	Bite Block, 20 per Box
------	------------------------



### Lavage Tubes

Ethox Lavage Tubes are Latex free and completely disposable. Tubes are 48" long with reference marks at 18", 22", and 20" from the distal end

2124	24fr Lavage Tube - 10 units per box (2lbs)
2128	28fr Lavage Tube - 10 units per box (2lbs)
2132	32fr Lavage Tube - 10 units per box (3lbs)
2136	36fr Lavage Tube - 10 units per box (3lbs)



**Ethox International, Inc.**

Tel: 1.800.521.1022

ethoxint.com

Fax: 716.842.2546